


## Map your proteins to OrthoMCL groups

This tool allows you to map proteins from a FASTA file to OrthoMCL groups. The tool has been implemented as a workflow in the VEuPathDB Galaxy workspace. To use this tool, you can follow these steps:

1. Click on the “Tools” item in the grey menu then select “Assign your proteins to groups” or click on the “Assign your proteins to groups” link from the right-hand side of the main OrthoMCL page.



2. The next page provides some instruction and a link to the VEuPathDB Galaxy server. Click on the VEuPathDB Galaxy Server link to access this service.  
**Note:** to use this service you will have to create a VEuPathDB account.



## Map your proteins to OrthoMCL groups

If you have a .fasta file with a set of proteins, you can map the proteins to OG5 or OG6r1 Groups.

If your .fasta file contains the proteins from a single proteome, you can additionally find paralog groups for proteins that do not map to OrthoMCL Core or Residual Groups.

First please log in as a VEuPathDB user.

Once logged in, go to the **VEuPathDB Galaxy server**:

1. Use the **Get Data** tool (in the left panel at the top) to upload your .fasta file
2. On the Galaxy home page, click **Workflow to map your proteins to OrthoMCL groups**.
3. In the workflow settings, indicate (1) your .fasta file, (2) the OG5 or OG6r1 Blast database, and (3) OG5 or OG6r1 groups.
4. Run the workflow!

**Note:** The Galaxy workflow run may take 24 hours or more to complete, depending on the size of the job queue

3. Once in Galaxy, you can import your protein FASTA file by clicking on the “Get Data” option in the left panel and selecting an import option, such as “Upload file from my computer”.

**globus Genomics** Analyze Data Workflow Shared Data Visualization Help User Using 485.9 GB

**Tools**

search tools

VEUPATHDB APPLICATIONS

- VEuPathDB Export Tools
- VEuPathDB OrthoMCL Tools
- VEuPathDB RNA-Seq Tools

DATA TRANSFER

- Globus Data Transfer
- Get Data**
- Upload File from your computer
- Get single dataset from BaseSpace
- Transfer dataset from BaseSpace by file ID
- Get Sample set Data from BaseSpace
- Transfer data from BaseSpace for a sample set

Collection Tools

- File Transfer Checksum

NGS APPLICATIONS

- NGS: QC and manipulation
- NGS: Assembly
- NGS: Mapping
- NGS: Mapping QC
- NGS: HLA Typing
- NGS: RNA Analysis
- NGS: miRNA
- NGS: DNase
- NGS: MTB

**Welcome to the VEuPathDB Galaxy Site**

A free, interactive, web-based platform for large-scale data analysis

**With VEuPathDB Galaxy you can:**

1. Start analyzing your data now with pre-configured workflows. All VEuPathDB genomes are pre-loaded.
2. Perform large-scale data analysis with no prior programming or bioinformatics experience.
3. Create custom workflows using an interactive workflow editor. [Learn how](#)
4. Export your results to VEuPathDB, so that you can explore your data with our tools, such as JBrowse and search strategies. See [this tutorial](#).
5. View your results on Galaxy or download results to your computer.
6. Keep data private, or share data with colleagues or the community.

To learn more about Galaxy, visit the [public Galaxy resources](#).

**Get started with VEuPathDB pre-configured workflows:**

**OrthoMCL**

This workflow uses BLASTP and the OrthoMCL algorithm to assign your set of proteins to OrthoMCL groups. Version OG6r1 is the latest set of groups (as of April 2020), but you can also select the previous set (OG5). [Explore this tutorial to learn more.](#)

- Workflow to map your proteins to OrthoMCL groups

**History**

search datasets

Unnamed history

(empty)

This history is empty. You can load your own data or get data from an external source

4. Your imported file will appear as a step in the right hand history panel. The color of the step indicates its status. The step is grey when it is in queue, yellow when in process and green when completed.
5. Once your protein FASTA file has been successfully uploaded into Galaxy, select the “Workflow to map your proteins to OrthoMCL groups” from the middle section.

**globus Genomics** Analyze Data Workflow Shared Data Visualization Help User Using 485.9 GB

**Tools**

search tools

VEUPATHDB APPLICATIONS

- VEuPathDB Export Tools
- VEuPathDB OrthoMCL Tools
- VEuPathDB RNA-Seq Tools

DATA TRANSFER

- Globus Data Transfer
- Get Data
- Collection Tools
- File Transfer Checksum

NGS APPLICATIONS

- NGS: QC and manipulation
- NGS: Assembly
- NGS: Mapping
- NGS: Mapping QC
- NGS: HLA Typing

**Get started with VEuPathDB pre-configured workflows:**

**OrthoMCL**

This workflow uses BLASTP and the OrthoMCL algorithm to assign your set of proteins to OrthoMCL groups. Version OG6r1 is the latest set of groups (as of April 2020), but you can also select the previous set (OG5). [Explore this tutorial to learn more.](#)

- Workflow to map your proteins to OrthoMCL groups

**RNA-seq**

Use the following workflows to analyze your FASTQ files. The workflows use FASTQ groomer and Trimmomatic for preparation of reads, FASTQC for sequencing statistics, and HISAT2 for mapping reads to a VEuPathDB reference genome. Choose the appropriate workflow based on your input data and your desired analysis. [Explore this tutorial](#) to learn about exporting your workflow results to VEuPathDB.

**History**

search datasets

Unnamed history

1 shown

5.9 MB

1: MyProteins

6. This will import the workflow into your workspace and allow you to run it on your protein file. In most cases you will not need to change any of the default parameters. Simply click on the “Run workflow” button at the bottom of the middle section. Note that there is a few second lag between clicking on the button and the workflow starting to run – please be patient.

### Running workflow “imported: Map your proteins to OrthoMCL groups”

[Expand All](#)[Collapse](#)

Assign your set of proteins to OrthoMCL groups. This workflow uses BLASTP and the OrthoMCL algorithm to (1) map proteins to OrthoMCL groups using BLAST similarity and (2) form paralog groups from proteins with no significant similarity to any OrthoMCL proteins. The workflow produces a file with the mapping from protein ID to group ID, along with similarity metrics. It also produces a file of paralog groups. The latter file is only valid if your input proteins all belong to a single proteome. This workflow might take 24 hours or more to run, depending on the size of the job queue.

[Step 1: OrthoMCL Clean FASTA file](#) (version 1.0.0)

1

**Proteome FASTA file:**

1: MyProteins 

**Maximum allowed number of input sequences**

100000 

**Action:**

Hide output 'output'.

[Step 2: NCBI BLAST+ makeblastdb](#) (version 0.3.0)

7

[Step 3: NCBI BLAST+ blastp](#) (version 0.3.0)

3

[Step 4: NCBI BLAST+ blastp](#) (version 0.3.0)

2

[Step 5: OrthoMCL Reformat Blast](#) (version 1.0.0)

5

[Step 6: OrthoMCL Reformat Blast](#) (version 1.0.0)

4

[Step 7: OrthoMCL Map Proteome to Groups](#) (version 1.0.0)

6

[Step 8: MCL](#) (version 14.137)

7

☐ Send results to a new history

[Run workflow](#)

7. Workflow steps will queue up in the right-hand panel. The entire workflow may take up to 24 hours or more to run depending on the size.

The screenshot shows the Galaxy Genomics interface. On the left is a sidebar with various tools categorized under 'Tools', 'Get Data', and 'NGS APPLICATIONS'. The main panel displays a green success message: 'Successfully ran workflow "Imported: Map your proteins to OrthoMCL groups". The following datasets have been added to the queue:'. Below this message is a list of 10 steps: 2: MyProteins.clean, 3: protein.BLAST.database.from.data.2, 4: blastp.MyProteins.clean.vs.orthomcl.v5.proteins.blast.db, 5: blastp.MyProteins.clean.vs.protein.BLAST.database.from.data.2, 6: reformat.blastp.MyProteins.clean.vs.cn\_sq.orthomcl.v5.proteins.blast.db\_sq.cn, 7: reformat.blastp.MyProteins.clean.vs.cn\_sq.protein.BLAST.database.from.data.2\_sq.cn, 9: paralogPairs.txt, 8: proteinsMappedToGroups.txt, and 10: MCS.on.data.9. On the right, the 'History' panel shows a list of datasets, including 'AssignToOrthoMCL' (11.91 MB) and 'MCS on data 9' (10 shown). The workflow steps are listed in a queue, with '1: MyProteins' at the bottom.

8. Once the workflow has completed, you can work with the output files. Click on the name of the step in your history in the right-hand panel to expand it. This provides additional options including the option to download file.

The screenshot shows the 'History' panel with the '10: MCS on data 9' step expanded. The expanded view shows the following information: '3 shown, 7 hidden', '78.32 MB', and a list of files: '10: MCS on data 9', '8: proteinsMappedToGroups.txt', '7,002 lines', 'format: tabular, database: ?', 'Scanning groups file /mnt/galaxyIndices2/genomes /miscellaneous /OrthoMCL\_Map\_Proteome\_to\_Groups.txt to build hash', 'Scanning similarity file /scratch/galaxy /files/051/dataset\_51041.dat', and 'Scanning self similarity file /scratch /galaxy/files/051/'. At the bottom, there is a 'Download' button and a list of files: 'mRNA1\_NF00000110-p1 OG5\_129304 ddisldhk', 'mRNA1\_NF0000020-p1 OG5\_197203 creil156', 'mRNA1\_NF0000030-p1 OG5\_159128 pchrifge', 'mRNA1\_NF0000070-p1 OG5\_179610 calbical', 'mRNA1\_NF0000090-p1 OG5\_126841 tpelest', and 'mRNA1\_NF0000110-p1 OG5\_187585 tniglens'. The '1: MyProteins' step is also visible at the bottom of the history panel.